MASKENOMICS® 101

Understanding ASTM Face Mask Performance Levels

Selecting the appropriate mask for a particular procedure is ensure your mask will provide the appropriate protection a critical component of your Personal Protective Equipment to minimize the spread of potentially infectious diseases. (PPE) protocol. Although masks may look similar, each mask has notable differences affecting the performance Medical face mask materials covered under ASTM and level of fluid resistance and particle filtration. Understanding the ASTM performance level of each face following performance levels based on the barrier mask can help make the selection process easier and

specifications are designated as one or more of the performance properties of the materials used:

ASTM F2100-11 STANDARDS

MASKENOMICS

performed.

Crosstex makes the selection process easy

filtration they provide... from minimum to

maximum performance. Just look to the

symbols displayed to find the appropriate

level of filtration for the procedure to be

with the MaskEnomics® rating system, designed to assist you in the decision

making process by categorizing all

Crosstex face masks by the level of

	ASTM Level 1	ASTM Level 2	ASTM Level 3
FLUID RESISTANCE, mmHg	80	120	160
BFE	≥ 95%	≥ 98%	≥ 98%
PFE, @ 0.1 micron	≥ 95%	≥ 98%	≥ 98%
DELTA P, mm H ₂ O/cm ²	< 4.0	< 5.0	< 5.0
FLAME SPREAD	Class 1	Class 1	Class 1

What the terms mean:

FLUID RESISTANCE represents the mask's resistance to penetration by synthetic blood under pressure (mmHg). It measures the ability of a mask's material construction to minimize fluids from traveling through the material and potentially coming into contact with the wearer. The higher the fluid resistance (filtration), the better the protection.

BFE (Bacterial Filtration Efficiency) represents the percentage of aerosol particulates filtered at a size of 3 microns. It is the measure of the efficiency of the mask in filtering bacteria passing through it.

PFE (sub-micron Particulate Filtration Efficiency) represents the percentage of submicron particulates filtered at 0.1 microns. PFE is the measure of the efficiency of the mask in filtering particles passing through it. The size of the particles filtered is critical.

DELTA P (Differential Pressure) represents the pressure drop across the mask or resistance to air flow in mmH₂O/cm². This determines breathing resistance – the higher the Delta P, the less the breathability, but the better the filtration.

FLAME SPREAD is a ranking derived by laboratory standard test methodology of a material's propensity to burn rapidly and spread flames.

SOURCE: ASTM Standard specification for performance of materials used in medical face masks – F2100-11 Standard.

Ordering Information

Isolite® Earloop Blue Pink	50/Box			/ D	
Pink	50/Box				
		GCLBL	Procedural Earloop	,	660011
	50/Box	GCLPK	Pink	50/Box	GCPPK
Yellow	50/Box	GCLYE	Yellow	50/Box	GCPYE
Surgical Molded			Pink Yellow Blue, Bulk Pack Surgical Tie-on	100/Box	GCPBL100
Blue	50/Box	GCBL	Surgical Tie-on		
Pink	50/Box 50/Box	GCPK	Blue	50/Box	GCS
Isofluid® Earloop w/ Se	cure Fit® Techno 50/Box			rloop w/ Secure Fit® Tecl	
Blue Lavender	50/Box 50/Box	GCIBLSF	White	50/Box	GCFCXSSI
Pink	50/Box 50/Box	GCILVSF	Ultra® Sensitive Foo	Free® Earloop	
		GCIPKSF	w/ Secure Fit® Techr	nology	
Isofluid® FogFree® Earl	•		White	40/Box	GCFCXSF5
w/ Secure Fit® Technolo	J,	CCICVECE	Ultra® Sensitive Foo	Free® Farloon	
Blue	40/Box	GCICXBSF	w/ Shield w/ Secure		
Isofluid® FogFree® Earlo	op w/ Shield		White	25/Box	GCPWSSF
w/ Secure Fit® Technolog	ıy				
Blue	25/Box	GCIPWBSF	-	ecure Fit® Technology	CCECVIIC
Isofluid® Earloop			Blue	50/Box	GCFCXUS
Blue	50/Box	GCIBL	Ultra® Sensitive Earl White Ultra® Sensitive Fog	loop	
Green	50/Box	GCIGR	White	50/Box	GCFCXS
Lavender	50/Box	GCILV	Ultra® Sensitive Fog	Eroe® Earloon	
Pink	50/Box	GCIPK	White	40/Box	GCFCXSF
Sapphire	50/Box	GCISA	1111111		- GCI CX51
Teal	50/Box	GCITE	_	Free® Earloop w/ Shield	
Topaz	50/Box	GCITZ	White	25/Box	GCPWS
Turquoise	50/Box	GCITQ	Ultra® Earloop		
Blue, Bulk Pack NEW!	100/Box	GCIBL100	Blue	40/Box	GCFCXU
Isofluid® FogFree® Earl	oop		Ultra® FogFree® Ear	loop	
Blue	40/Box	GCICXB	Blue	40/Box	GCFCX
Sapphire	40/Box	GCICXS	Шене в Бенебие в Вене	la au/ Chialal	
Topaz	40/Box	GCICXZ	Ultra® FogFree® Ear Blue		GCPW
Turquoise	40/Box	GCICXT		25/Box	GCPW
Isofluid® FogFree® Earl	oon w/ Shield		6)	articulate Respirator	
Blue	25/Box	GCIPWB	Blue/white stripe	28/Box	GPRN95
Sapphire	25/Box	GCIPWS	Face Shield		
Topaz	25/Box 25/Box	GCIPWZ	Clear	24/Box	GCSS
Turquoise	25/Box 25/Box	GCIPWT	Clear, Bulk Pack	72/Box	GCSSB
<u> </u>	25/00%	GCII VV I			
Isofluid® Plus Earloop	50/D	CDITICDI	Patient Safety Mask White	25/Box	GCPAT
Blue	50/Box	GPLUSBL			GCFAI
Kaleidoscope	50/Box	GPLUSKA	See-Clear™ Eyeglas		
Pink	50/Box	GPLUSPK	Wipes	120/Box	SEEC
Isofluid® Plus Pure w/ S	ecure Fit® Techn	ology			
White NEW!	50/Box	GPLUSWHSF	= NOT MADE WITH NAT	TIRAI RIIRRERI ATEX	
Procedural Earloop w/	Secure Fit® Tech	nology	- NOT MADE WITH NAT	ONAL NOBBEN LATEA	



10 Ranick Road Hauppauge, NY 11788-4209 USA

Tel: 631-582-6777 | 888-276-7783 Fax: 631-582-1726 | E-mail: crosstex@crosstex.com

California • Colorado • Georgia • New York • Pennsylvania • Argentina • Japan • Singapore • The Netherlands

Crosstex can make no warranties or representations, either expressed or implied, that these products will fully protect the user from exposure to blood or bodily fluids or risk of contracting infectious diseases. OSHA requires the employer to evaluate the anticipated exposure and select the appropriate protective masks to prevent contamination of skin, eyes and respiratory passages.

Lavender

Lavender

Procedural Earloop

Pink



GCPBLSF

GCPLVSF

GCPPKSF

GCPBL

GCPLV

















50/Box

50/Box

50/Box

50/Box

WHY CROSSTEX®?

Enhanced Safety and Performance!

Crosstex offers a wide variety of mask designs, fit and filtration to match the protection needs for each procedure or risk level. The Crosstex MaskEnomics® filtration guide makes it easy to select the level of filtration required, including ASTM Level 1, 2 and 3 and N95 particulate respiration. From minimum performance masks for dry, short procedures, to our maximum filtration mask, Isolator Plus, Crosstex has the right mask for the task!





SECÜRE FIT. Technology

TOP MASK 3 YEARS IN A ROW Experience Secure Fit® for vourself. **Better Fit.**

Better Filtration. Better Protection.

Although the N95 closely resembles a surgical mask, it is a

respirator and must be used in accordance with all

OSHA regulations regarding respiratory protection.

3X More Protection*

- → Aluminum strips both above the nose and below the chin for a custom fit. Conforms to any size or shape face.
- → Patent pending design significantly reduces gapping at top, bottom and sides.

*Study on File



FEATURE GUIDE TO CROSSTEX EARLOOP/TIE-ON FACE MASKS A B C D E F G H I J K L M tra® Sensitive FogFree® Earloop w/ Secure Fit® Technolog Ultra® Sensitive Earloop Ultra® Sensitive FogFree® Earloop Ultra® Sensitive FogFree® Earloop w/ Shield Itra® FogFree® Earloon Jltra® FogFree® Earloop w/ Shield

B Spunbond outer layer material for naximum protection against fluids. Facial tissue outer layer is ultrasoft

maximum protection.

A Extra-long, enclosed aluminum

nose piece forms a strong seal for

D Vapor barrier uniquely placed on outside layer to prevent fogging.

E FogFree / No-Fog strip blocks and orbs moisture; forms a strong seal to prevent fogging; cushions nose ridge.

..................

OUTSIDE

F Secure Loop earloops are soft, strong, t comfortabley, and attached to outside of mask to eliminate irritation.


G Soft-Bond ultrasonics ensures strongest construction, eliminating holes and defects.

H Highest quality filter media for best ASTM rating, filtration value and breath-

inks and chemicals, which can cause

Cellulose inner layer for sensitive skin (white, extra soft, hypoallergenic will not lint, tear or shred; void of dyes,

J Tissue inner layer is made of white, oft facial tissue.



K Spunbond inner layer made of white, luid-resistant material.

Attached wrap-around face shield optically clear and distortion-free.

Secure Fit® aluminum strip under chin significantly reduces gapping at the top. sides and bottom of the mask.

All Crosstex face masks are not made with natural rubber latex.

ASTM LEVEL 3

Ideal for procedures where heavy to moderate amounts of fluid, spray and/or aerosols are produced.

MAXIMUM FILTRATION

Meets EN 149: 2001 FFP2 NR

Meets EN14683 Rating - Type IIR Standard.

Indicated for use when treating patients with

airborne diseases such as TB or influenza.



Ultra® Sensitive Earloop w/ Secure Fit® Technology

- BFE = 99.9% @ 3 microns.
- for custom fit. Fluid resistant outer laver.

- PFE = 99.8% @ 0.1 micron.
- Patent pending design with aluminum strip under chin

Isofluid® Earloop w/ Secure Fit® Technology

Isolator Plus N95 Particulate Respirator

Meets guidelines for TB exposure control

Resistance to penetration by synthetic blood:

NIOSH approved N95 Respirator.

White hypoallergenic inner cellulose

· Malleable nose and chin pieces.

• Inner/outer fluid resistant layers.

• Extra soft – will not lint, tear or shred. Free of chemicals, inks and dyes. • Extra-long aluminum nose piece.



ÜŠÄ 🕾 C E



Ultra® Sensitive FogFree® Earloop w/ Secure Fit® Technology

- BFE = 99.9% @ 3 microns. PFE = 99.8% @ 0.1 micron.
- Provides all the features of Ultra® Sensitive with Secure Fit® technology with FogFree® strip üšä 🕪 (E

Ultra® Sensitive FogFree® Earloop w/Shield w/ Secure Fit®Technology

- BFE = 99.9% @ 3 microns.
- PFF = 99 8% @ 0.1 micron Provides all the features of Ultra®Sensitive FogFree® with Secure Fit® technology with attached optically clear.

distortion-free shield. () USA (E

Ultra® Earloop w/ Secure Fit®Technology • BFE = 99.9% @ 3 microns.

- PFE = 99.8% @ 0.1 micron.
- Patent pending design with aluminum

Isofluid® FogFree® Earloop

w/ Shield w/ Secure Fit®

Technology

- strip under chin for tight, custom fit. Fluid resistant white inner layer. Extra-long aluminum nose piece
- Fluid resistant outer layer. 👸 🐯 🍖 🕻 🧲

Ultra® Sensitive Earloop

- BFE = 99.9% @ 3 microns.
- PFE = 99.8% @ 0.1 micron.
- Fluid resistant outer layer. White hypoallergenic inner cellulose layer
- Extra soft will not lint, tear or shred.
- Free of chemicals, inks and dyes. • Extra-long aluminum nose piece.

SAM CE Also available: Ultra® Sensitive FogFree® Earloop Ultra® Sensitive FogFree® Earloop with

• BFE = 99.9% @ 3 microns. • PFF = 99 8% @ 0.1 micron

Fluid resistant outer laver.

Also available:

Ultra® FogFree® Earloop

Ultra® Earloop

 Fluid resistant white inner laver. · Extra-long aluminum nose piece.

Ultra® FogFree® Earloop with Shield



ward Winning Crosstex Masl

2010

w/ Shield

· Unique FogFree® with optically clear, distortion free wrap around shield

Ultra® FogFree® Earloop

• BFE = >99.9% @ 3 microns

üsä 🌬 (E

ASTM LEVEL 2

ASTM LEVEL 1

Ideal for procedures where moderate to light amounts of fluid, spray and/or aerosols are produced.

Ideal for procedures where low amounts of

fluid, spray and/or aerosols are produced.

Ideal as a comfortable substitute for earloop

procedures that do not produce fluid, spray

face masks, this mask is a simple physical barrier

ideal for exams and visitations or for dry, short

Meets EN14683 Rating - Type II Standard.

LOW PERFORMANCE

or aerosols.

Meets EN14683 Rating – Type IIR Standard.



Procedural Earloop w/ Secure Fit® Technology

- BFE = 99.9% @ 3 microns. PFE = 99.7% @ 0.1 micron.
- Patent pending design with aluminum strip under chin for custom fit.

Patent pending design with aluminum strip under chin for custom fit.

- Fluid resistant outer laver.
- Fluid resistant white inner layer
- Extra-long aluminum nose piece.

• BFE = 98.4% @ 3 microns.

PFE ≥ 95% @ 0.1 micron.

Fluid resistant outer layer.

Extra-long aluminum nose piece.

White tissue inner layer.

• BFE = 98.7% @ 4.2 microns.

· Flared edge prevents irritation.

Splash resistant barrier.

· Soft, flexible nose piece.



JUSA 🕪 (E

USA 🏁 🤄 (E

ÜSÄ 🖗 C E

Procedural Earloop

- BFE = 99.9% @ 3 microns. • PFE = 99.7% @ 0.1 micron.
- Fluid resistant outer layer. Fluid resistant white inner layer.
- Extra-long aluminum nose piece.

Isofluid® Earloop

• BFE = 98.4% @ 3 microns.

• PFE ≥ 95% @ 0.1 micron.

· White tissue inner layer.

Fluid resistant outer layer.

Extra-long aluminum nose piece



USA 🗪 (E

Surgical Tie-on BFE = 98% @ 1 micron.

- PFE = 99.2% @ 0.1 micron.
- Fluid resistant outer layer. Fluid resistant white inner layer. Extra-large panel for surgical procedures.

Isofluid® FogFree® Earloop

w/ Secure Fit® Technology

Provides all features of Secure Fit[®]

BFF = 98 4% @ 3 microns

• PFE ≥ 95% @ 0.1 micron.

Isofluid with addition of

Extra-long aluminum nose piece.

Isofluid® Plus Earloop

- BFE = 98.4% @ 3 microns. • PFE ≥ 95% @ 0.1 micron. Provides all features of Secure Fit® Isofluid FogFree® with attached optically clear ③ üšä⋙ 🚱 🤇 € distortion-free



• BFE = 99.8% @ 3 microns.

- PFE = 99.01% @ 0.1 micron.
- Fluid resistant outer layer. · White hypoallergenic inner cellulose la · Extra-long aluminum nose piece. • Extra soft – will not lint, tear or shred.





Isofluid® Plus Pure w/Secure Fit® Technology

- BFE = 99.8% @ 3 microns. PFE = 99.01% @ 0.1 micron. Fluid resistant outer layer.
- · White hypoallergenic inner cellulose layer. · Extra-long aluminum nose piece. • Extra soft – will not lint, tear or shred.



USA 🗪 CE



The Dental Advisor: **▶** Secure Fit® 2012

2013

2014

2008 2011

2009 2012

ADVISOR

▶ Patient Safety Mask with Shield ++++

The Dental Advisor:

▶ Innovative **Products**

Patient Protection

MINIMUM PERFORMANCE

Ideal as a simple physical barrier for exams, isolation and visitations or for dry, short procedures that do not produce fluid, spray or aerosols.



Surgical Molded

Isolite® Earloop • BFE = 99.6% @ 3 microns.

 Outer facial tissue laver. White facial tissue inner laver Extra-long aluminum nose piece.

№ C €

with no pressure on temples; vented for increased air flow.

 Sonically welded elastic headband for added strength. Anti-fog treatment on inside and outside of shield.

USA (E

Additional Face Protection Products

The mask performance levels were performed for Crosstex by Nelson Laboratories in accordance with ASTM Standard F2100-11 and represent the most current performance test

Isofluid FogFree® Earloop with Shield

Also available:

Isofluid FogFree® Earloop

See-Clear™ **Eyeglass Wipes**

- 4" x 5.5" wipes are anti-fog,
- · Will not scratch, dries quickly.


Patient Safety Mask with Shield

- U.S. Patent No.'s 6.185,740 and 6.523,170.
- Covers patient's eyes and nose. Protects fragile eye and mucosal tissue from flying particles, spray and
- Makes patients feel safe, but not "in the dark" with attached clear-vision shield.
- Form-fitting profile, flexible materials provide maximum access to mouth. • White spunbond inner and outer layers.
- Fluid resistant.









- · Protects mask and face from direct splatter; may prolong mask life.

• Effective cleaning for eyewear and face shields.



anti-static.

