



PERSONAL PROTECTION



PERSONAL PROTECTION

### LATEX EXAM GLOVES-AMBIDEXTROUS

#### Exam - Lightly Powdered

- Smooth, natural/white.
- Powder: USP modified corn starch inside, USP Calcium carbonate outside.



100/Box

Description	Ref. #	Description	Ref. #
<input type="checkbox"/> Extra Small	19003	<input type="checkbox"/> Medium	19001
<input type="checkbox"/> Small	19000	<input type="checkbox"/> Large	19002

#### Omega - Lightly Powdered

- Textured, natural/white.
- Powder: USP modified corn starch inside, USP calcium carbonate outside.



100/Box

Description	Ref. #	Description	Ref. #
<input type="checkbox"/> Extra Small	19080	<input type="checkbox"/> Medium	19083
<input type="checkbox"/> Small	19081	<input type="checkbox"/> Medium/Large	19084
<input type="checkbox"/> Small/Medium	19082	<input type="checkbox"/> Large	19085

#### Classic® - Lightly Powdered

- Smooth, bisque finish.
- Allotex – latex with extremely low antigenic protein.
- Low modulus – enhanced comfort.
- Powder: USP modified cornstarch.



100/Box

Description	Ref. #	Description	Ref. #
<input type="checkbox"/> Small	MSC	<input type="checkbox"/> Large	MLC
<input type="checkbox"/> Medium	MMC		

### LATEX EXAM GLOVES-AMBIDEXTROUS

#### Ultra Plus® - Powder-Free

- Smooth, bisque-finish, no pigment.
- Polymer lined, non-chlorinated, double washed inside & out.



100/Box

Description	Ref. #	Description	Ref. #
<input type="checkbox"/> Small	MSCP	<input type="checkbox"/> Large	MLCP
<input type="checkbox"/> Medium	MMCP		

#### Friction Grip - Powder-Free

- Textured, natural/white.
- Double polymer coated for donning.



100/Box

Description	Ref. #	Description	Ref. #
<input type="checkbox"/> Extra Small	19043	<input type="checkbox"/> Medium	19041
<input type="checkbox"/> Small	19040	<input type="checkbox"/> Large	19042

For the full line of Crosstex products visit [crosstex.com](http://crosstex.com).  
To order, please contact your authorized Crosstex dealer.

**CROSSTEX** | PROTECTS®  
A CANTEL MEDICAL COMPANY

Manufactured by or for Crosstex  
10 Ranick Road, Hauppauge, NY 11788-4209 USA  
Tel: 631-582-6777 | 888-276-7783 | Fax: 631-582-1726  
E-mail: [crosstex@crosstex.com](mailto:crosstex@crosstex.com) | [crosstex.com](http://crosstex.com)

California • Colorado • Georgia • New York • Pennsylvania • Argentina • Japan • Singapore • The Netherlands  
Proud member of the following associations:



## YOUR SAFETY IS IN YOUR HANDS!



Crosstex manufactures a majority of their products in the U.S.A.

Click-N-Save at [crosstex.com](http://crosstex.com)



888-276-7783  
[crosstex.com](http://crosstex.com)

**CROSSTEX** | PROTECTS®  
A CANTEL MEDICAL COMPANY

© 2014 Crosstex - 2/11/14 - CTX1143



## PERSONAL PROTECTION

Gloves are the most important personal barrier to protect against blood-borne pathogen exposure. Every office should have a selection of gloves to meet the needs of each procedure and staff member.

### FITTED LATEX EXAM GLOVES-NON-STERILE

#### Left/Right Fitted - Lightly Powdered

- Textured, no color pigment.
- Powder: USP modified corn starch.



50 Pair/Box

Description	Ref. #	Description	Ref. #
<input type="checkbox"/> 6.5	19031	<input type="checkbox"/> 8	19034
<input type="checkbox"/> 7	19032	<input type="checkbox"/> 8.5	19035
<input type="checkbox"/> 7.5	19033	<input type="checkbox"/> 9	19036

#### Ultimate Left/Right Powder-Free

- Textured, natural/white.
- Double polymer coated for donning.



50 Pair/Box

Description	Ref. #	Description	Ref. #
<input type="checkbox"/> 6.5	19071	<input type="checkbox"/> 8	19074
<input type="checkbox"/> 7	19072	<input type="checkbox"/> 8.5	19075
<input type="checkbox"/> 7.5	19073	<input type="checkbox"/> 9	19076

### FITTED NON-LATEX GLOVES

#### Nitrile Left/Right Fitted Powder-Free

- Textured, blue 100% synthetic nitrile.
- Double polymer coated for donning.



50 Pair/Box

Description	Ref. #	Description	Ref. #
<input type="checkbox"/> 6.5	19093	<input type="checkbox"/> 8	19096
<input type="checkbox"/> 7	19094	<input type="checkbox"/> 8.5	19097
<input type="checkbox"/> 7.5	19095		

#### Copolymer Plastic Powder-Free

- Smooth, 100% clear synthetic polymer.
- Overglove, not for patient treatment.
- Dispensed in individual sheets.



100/Box

Description	Ref. #	Description	Ref. #
<input type="checkbox"/> Small	MSCC	<input type="checkbox"/> Large	MLCC
<input type="checkbox"/> Medium	MMCC		

**NEW!** *Ambi-ance™*

#### Ambidextrous Nitrile Glove Powder-Free

- Soft stretch nitrile.
- Puncture and tear resistant barrier.
- Textured fingertips.



200/Box

Description	Ref. #	Description	Ref. #
<input type="checkbox"/> Extra Small	19010	<input type="checkbox"/> Medium	19012
<input type="checkbox"/> Small	19011	<input type="checkbox"/> Large	19013

All Crosstex gloves meet FDA, ASTM and ISO physical and performance regulations and specifications for medical grade gloves, including:

- Physical dimensions, thickness, tensile strength and modulus.
- Major defects (pinholes, tears), minor defects (visual).

GLOVE	AMBIDEXTRIOUS LATEX EXAM GLOVES					FITTED LATEX GLOVES		NON-LATEX GLOVES (CONTAINS NO LATEX PROTEIN)		
	EXAM	OMEGA	CLASSIC	ULTRA PLUS	FRICTION GRIP	FITTED LATEX	ULTIMATE FITTED	NITRILE FITTED	AMBIDEXTRIOUS NITRILE	COPOLYMER PLASTIC
POWDERED	•	•	•			•				
POWDER-FREE				•	•		•	•	•	•
MEETS ASTM D5712 MAX. TOTAL PROTEIN OF < 200 µg/dm	•	•	•	•	•	•	•	N/A	N/A	N/A
POLYMER-COATED			•	•	Inside and Out	Inside and Out	Inside and Out	Inside and Out		
NON-CHLORINATED	•	•	•	•	•	•	•	•	N/A	N/A
TEXTURE	Smooth	Textured	Smooth	Smooth	Textured	Textured	Textured	Textured	Textured Fingers	Smooth
COLOR	Natural White	Natural White	Natural White	Natural White	Natural White	Natural White	Natural White	Blue	Blue	Clear
BEADED CUFF	•	•	•	•	•	•	•	•	•	N/A

= NOT MADE WITH NATURAL RUBBER LATEX