

Sterilization Compliance Calendar

Keep track of sterilization procedures to protect your practice & your patients!

2018

For customer test results visit:
www.CrosstexBMS.com

Or call Customer Care at:
1-800-819-3336

How often should I perform biological monitoring?

CDC and AAMI guidelines recommend biological testing should be performed on a weekly basis.

Please check state recommendations for additional information.

Visit: cdc.gov or aami.org

In addition to weekly monitoring, practices should perform biological monitoring:

- ✓ Whenever a new type of packaging material or tray is used
- ✓ After training new sterilization personnel
- ✓ After any sterilizer repair or maintenance
- ✓ After any change in the sterilizer loading process

HELPFUL REMINDERS



Never run the BI on the first load of the day



Always run the BI on a full load if using a table top sterilizer



Sterilant penetration best monitored by placing BIs in the hardest to reach areas of sterilizer

KEY



Place your initials over the date each week that BI testing was performed.



EVERY LOAD. SURE-CHECK[®] Sterilization Pouches Indicators will turn from pink to brown indicating Time/Temp/Steam sterilization requirements have been met.



EVERY DAY. STEAMPlus[™] Type 5 Integrators Use a check mark to indicate daily testing with a Type 5 Integrator. Moving-front technology must reach into SAFE zone to pass.



EVERY WEEK. CONFIRM[®] Biological Monitoring System Indicate Pass (P), Fail (F) or Invalid (I).

JANUARY

SUN	MON	TUE	WED	THU	FRI	SAT	WEEKLY RESULTS (P/F/I)
	1	2	3	4	5	6	
7	8	9	10	11	12	13	
14	15	16	17	18	19	20	
21	22	23	24	25	26	27	
28	29	30	31				

FEBRUARY

SUN	MON	TUE	WED	THU	FRI	SAT	WEEKLY RESULTS (P/F/I)
				1	2	3	
4	5	6	7	8	9	10	
11	12	13	14	15	16	17	
18	19	20	21	22	23	24	
25	26	27	28				

MARCH

SUN	MON	TUE	WED	THU	FRI	SAT	WEEKLY RESULTS (P/F/I)
				1	2	3	
4	5	6	7	8	9	10	
11	12	13	14	15	16	17	
18	19	20	21	22	23	24	
25	26	27	28	29	30		

APRIL

SUN	MON	TUE	WED	THU	FRI	SAT	WEEKLY RESULTS (P/F/I)
1	2	3	4	5	6	7	
8	9	10	11	12	13	14	
15	16	17	18	19	20	21	
22	23	24	25	26	27	28	
29	30						

MAY

SUN	MON	TUE	WED	THU	FRI	SAT	WEEKLY RESULTS (P/F/I)
		1	2	3	4	5	
6	7	8	9	10	11	12	
13	14	15	16	17	18	19	
20	21	22	23	24	25	26	
27	28	29	30	31			

JUNE

SUN	MON	TUE	WED	THU	FRI	SAT	WEEKLY RESULTS (P/F/I)
					1	2	
3	4	5	6	7	8	9	
10	11	12	13	14	15	16	
17	18	19	20	21	22	23	
24	25	26	27	28	29	30	

JULY

SUN	MON	TUE	WED	THU	FRI	SAT	WEEKLY RESULTS (P/F/I)
1	2	3	4	5	6	7	
8	9	10	11	12	13	14	
15	16	17	18	19	20	21	
22	23	24	25	26	27	28	
29	30	31					

AUGUST

SUN	MON	TUE	WED	THU	FRI	SAT	WEEKLY RESULTS (P/F/I)
			1	2	3	4	
5	6	7	8	9	10	11	
12	13	14	15	16	17	18	
19	20	21	22	23	24	25	
26	27	28	29	30	31		

SEPTEMBER

SUN	MON	TUE	WED	THU	FRI	SAT	WEEKLY RESULTS (P/F/I)
						1	
2	3	4	5	6	7	8	
9	10	11	12	13	14	15	
16	17	18	19	20	21	22	
23	24	25	26	27	28	29	
30							

OCTOBER

SUN	MON	TUE	WED	THU	FRI	SAT	WEEKLY RESULTS (P/F/I)
	1	2	3	4	5	6	
7	8	9	10	11	12	13	
14	15	16	17	18	19	20	
21	22	23	24	25	26	27	
28	29	30	31				

NOVEMBER

SUN	MON	TUE	WED	THU	FRI	SAT	WEEKLY RESULTS (P/F/I)
				1	2	3	
4	5	6	7	8	9	10	
11	12	13	14	15	16	17	
18	19	20	21	22	23	24	
25	26	27	28	29	30		

DECEMBER

SUN	MON	TUE	WED	THU	FRI	SAT	WEEKLY RESULTS (P/F/I)
						1	
2	3	4	5	6	7	8	
9	10	11	12	13	14	15	
16	17	18	19	20	21	22	
23	24	25	26	27	28	29	
30	31						