

2020

www.CrosstexBMS.com

Contact us for test status:
1-800-819-3336

Because you can't see sterile.™

Biological Monitoring Quarterly Compliance Calendar

How often should I perform biological monitoring (BI) (spore testing)?

The CDC recommends testing at least weekly with a biological monitor for each sterilizer to ensure proper functioning. Users should follow the manufacturer's directions concerning the appropriate placement of the BI in the sterilizer. In addition to conducting routine biological monitoring, equipment users should perform biological monitoring:

- Whenever a new type of packaging material or tray is used.
- After training new sterilization personnel.
- After any sterilizer repair or maintenance.
- After any change in the sterilizer loading procedures.

Circle the date each week that BI testing was performed

JANUARY

S	M	T	W	T	F	S
			1	2	3	4
5	6	7	8	9	10	11
12	13	14	15	16	17	18
19	20	21	22	23	24	25
26	27	28	29	30	31	

FEBRUARY

S	M	T	W	T	F	S
						1
2	3	4	5	6	7	8
9	10	11	12	13	14	15
16	17	18	19	20	21	22
23	24	25	26	27	28	29

MARCH

S	M	T	W	T	F	S
1	2	3	4	5	6	7
8	9	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28
29	30	31				

TEST DATE	INITIALS OF PERSON WHO RAN TEST	MAIL DATE	INITIALS OF PERSON WHO MAILED TEST

Have any questions?
Call Crosstex Customer Care at:
1-800-819-3336

FOR STERILITY ASSURANCE:
EVERY LOAD. Sure-Check™ Sterilization Pouches
EVERY DAY. STEAMPlus™ Type 5 Integrators
EVERY WEEK. ConFirm™ Biological Monitoring Services



All product names are trademarks of Crosstex International, Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted.



2020

www.CrosstexBMS.com

Contact us for test status:
1-800-819-3336

Because you can't see sterile.™

Biological Monitoring Quarterly Compliance Calendar

How often should I perform biological monitoring (BI) (spore testing)?

The CDC recommends testing at **least weekly** with a biological monitor for each sterilizer to ensure proper functioning. Users should follow the manufacturer's directions concerning the appropriate placement of the BI in the sterilizer. In addition to conducting routine biological monitoring, equipment users should perform biological monitoring:

- Whenever a new type of packaging material or tray is used.
- After training new sterilization personnel.
- After any sterilizer repair or maintenance.
- After any change in the sterilizer loading procedures.

Circle the date each week that BI testing was performed

APRIL

S	M	T	W	T	F	S
			1	2	3	4
5	6	7	8	9	10	11
12	13	14	15	16	17	18
19	20	21	22	23	24	25
26	27	28	29	30		

MAY

S	M	T	W	T	F	S
					1	2
3	4	5	6	7	8	9
10	11	12	13	14	15	16
17	18	19	20	21	22	23
24	25	26	27	28	29	30
31						

JUNE

S	M	T	W	T	F	S
	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30				

TEST DATE	INITIALS OF PERSON WHO RAN TEST	MAIL DATE	INITIALS OF PERSON WHO MAILED TEST

Have any questions?

Call Crosstex Customer Care at:

1-800-819-3336

FOR STERILITY ASSURANCE:

EVERY LOAD. Sure-Check™ Sterilization Pouches

EVERY DAY. STEAMPlus™ Type 5 Integrators

EVERY WEEK. ConFirm™ Biological Monitoring Services



All product names are trademarks of Crosstex International, Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted.

© 2019 Crosstex International, Inc. DLIT01341 Rev A 1219

(888) 276-7783

crosstex.com



2020

www.CrosstexBMS.com

Contact us for test status:
1-800-819-3336

Because you can't see sterile.™

Biological Monitoring Quarterly Compliance Calendar

How often should I perform biological monitoring (BI) (spore testing)?

The CDC recommends testing **at least weekly** with a biological monitor for each sterilizer to ensure proper functioning. Users should follow the manufacturer's directions concerning the appropriate placement of the BI in the sterilizer. In addition to conducting routine biological monitoring, equipment users should perform biological monitoring:

- Whenever a new type of packaging material or tray is used.
- After training new sterilization personnel.
- After any sterilizer repair or maintenance.
- After any change in the sterilizer loading procedures.

Circle the date each week that BI testing was performed

JULY

S	M	T	W	T	F	S
			1	2	3	4
5	6	7	8	9	10	11
12	13	14	15	16	17	18
19	20	21	22	23	24	25
26	27	28	29	30	31	

AUGUST

S	M	T	W	T	F	S
						1
2	3	4	5	6	7	8
9	10	11	12	13	14	15
16	17	18	19	20	21	22
23	24	25	26	27	28	29
30	31					

SEPTEMBER

S	M	T	W	T	F	S
		1	2	3	4	5
6	7	8	9	10	11	12
13	14	15	16	17	18	19
20	21	22	23	24	25	26
27	28	29	30			

TEST DATE	INITIALS OF PERSON WHO RAN TEST	MAIL DATE	INITIALS OF PERSON WHO MAILED TEST

Have any questions?
Call Crosstex Customer Care at:
1-800-819-3336

FOR STERILITY ASSURANCE:
EVERY LOAD. Sure-Check™ Sterilization Pouches
EVERY DAY. STEAMPlus™ Type 5 Integrators
EVERY WEEK. ConFirm™ Biological Monitoring Services



All product names are trademarks of Crosstex International, Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted.
© 2019 Crosstex International, Inc. DLIT01341 Rev A 1219

(888) 276-7783 | crosstex.com



2020

www.CrosstexBMS.com

Contact us for test status:
1-800-819-3336

Because you can't see sterile.™

Biological Monitoring Quarterly Compliance Calendar

How often should I perform biological monitoring (BI) (spore testing)?

The CDC recommends testing **at least weekly** with a biological monitor for each sterilizer to ensure proper functioning. Users should follow the manufacturer's directions concerning the appropriate placement of the BI in the sterilizer. In addition to conducting routine biological monitoring, equipment users should perform biological monitoring:

- Whenever a new type of packaging material or tray is used.
- After training new sterilization personnel.
- After any sterilizer repair or maintenance.
- After any change in the sterilizer loading procedures.

Circle the date each week that BI testing was performed

OCTOBER

S M T W T F S

					1	2	3
4	5	6	7	8	9	10	
11	12	13	14	15	16	17	
18	19	20	21	22	23	24	
25	26	27	28	29	30	31	

NOVEMBER

S M T W T F S

1	2	3	4	5	6	7	
8	9	10	11	12	13	14	
15	16	17	18	19	20	21	
22	23	24	25	26	27	28	
29	30						

DECEMBER

S M T W T F S

		1	2	3	4	5	
6	7	8	9	10	11	12	
13	14	15	16	17	18	19	
20	21	22	23	24	25	26	
27	28	29	30	31			

TEST DATE	INITIALS OF PERSON WHO RAN TEST	MAIL DATE	INITIALS OF PERSON WHO MAILED TEST

Have any questions?
Call Crosstex Customer Care at:
1-800-819-3336

FOR STERILITY ASSURANCE:
EVERY LOAD. Sure-Check™ Sterilization Pouches
EVERY DAY. STEAMPlus™ Type 5 Integrators
EVERY WEEK. ConFirm™ Biological Monitoring Services



All product names are trademarks of Crosstex International, Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted.

© 2019 Crosstex International, Inc. DLIT01341 Rev A 1219

