

# 2020

For customer test results visit:  
www.CrosstexBMS.com

Or call Customer Care at:  
1-800-819-3336

Because you can't see sterile.™

# Sterilization Compliance Calendar

Keep track of sterilization procedures to protect your practice & your patients!

**How often should I perform biological monitoring?**

CDC and AAMI guidelines recommend biological testing should be performed on a weekly basis.

Please check state recommendations for additional information.

Visit: [cdc.gov](http://cdc.gov) or [aami.org](http://aami.org)

**In addition to weekly monitoring, practices should perform biological monitoring:**

- ✓ Whenever a new type of packaging material or tray is used
- ✓ After training new sterilization personnel
- ✓ After any sterilizer repair or maintenance
- ✓ After any change in the sterilizer loading process

**HELPFUL REMINDERS**



Never run the BI on the first load of the day



Always run the BI on a full load if using a table top sterilizer



Sterilant penetration best monitored by placing BIs in the hardest to reach areas of sterilizer

**KEY**



Place your initials over the date each week that BI testing was performed.



**EVERY LOAD. Sure-Check™ Sterilization Pouches**  
Indicators will turn from pink to brown indicating Time/Temp/Steam sterilization requirements have been met.



**EVERY DAY. STEAMPlus™ Type 5 Integrators**  
Use a check mark to indicate daily testing with a Type 5 Integrator. Moving-front technology must reach into SAFE zone to pass.



**EVERY WEEK. ConFirm™ Biological Monitoring System**  
Indicate Pass (P), Fail (F) or Invalid (I).

**JANUARY**

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)
			1	2	3	4	
5	6	7	8	9	10	11	
12	13	14	15	16	17	18	
19	20	21	22	23	24	25	
26	27	28	29	30	31		

**FEBRUARY**

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)
						1	
2	3	4	5	6	7	8	
9	10	11	12	13	14	15	
16	17	18	19	20	21	22	
23	24	25	26	27	28	29	

**MARCH**

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)
1	2	3	4	5	6	7	
8	9	10	11	12	13	14	
15	16	17	18	19	20	21	
22	23	24	25	26	27	28	
29	30	31					

**APRIL**

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)
			1	2	3	4	
5	6	7	8	9	10	11	
12	13	14	15	16	17	18	
19	20	21	22	23	24	25	
26	27	28	29	30			

**MAY**

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)
					1	2	
3	4	5	6	7	8	9	
10	11	12	13	14	15	16	
17	18	19	20	21	22	23	
24	25	26	27	28	29	30	
31							

**JUNE**

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)	
		1	2	3	4	5	6	
7	8	9	10	11	12	13		
14	15	16	17	18	19	20		
21	22	23	24	25	26	27		
28	29	30						

**JULY**

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)
			1	2	3	4	
5	6	7	8	9	10	11	
12	13	14	15	16	17	18	
19	20	21	22	23	24	25	
26	27	28	29	30	31		

**AUGUST**

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)
						1	
2	3	4	5	6	7	8	
9	10	11	12	13	14	15	
16	17	18	19	20	21	22	
23	24	25	26	27	28	29	
30	31						

**SEPTEMBER**

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)
		1	2	3	4	5	
6	7	8	9	10	11	12	
13	14	15	16	17	18	19	
20	21	22	23	24	25	26	
27	28	29	30				

**OCTOBER**

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)
				1	2	3	
4	5	6	7	8	9	10	
11	12	13	14	15	16	17	
18	19	20	21	22	23	24	
25	26	27	28	29	30	31	

**NOVEMBER**

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)
1	2	3	4	5	6	7	
8	9	10	11	12	13	14	
15	16	17	18	19	20	21	
22	23	24	25	26	27	28	
29	30						

**DECEMBER**

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)
		1	2	3	4	5	
6	7	8	9	10	11	12	
13	14	15	16	17	18	19	
20	21	22	23	24	25	26	
27	28	29	30	31			



All product names are trademarks of Crosstex International, Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted.

© 2019 Crosstex International, Inc. DLIT01340 Rev A 1219

(888) 276-7783 [crosstex.com](http://crosstex.com)

