

2026

For customer test results visit:

CrosstexBMS.com

Or call Customer Care at:

1-800-722-1529

Sterilization Compliance Calendar

Keep track of sterilization monitoring procedures in your practice!

How often should I perform biological monitoring?

CDC and AAMI guidelines recommend biological testing should be performed at least weekly for each sterilizer to ensure proper functioning.

Please check individual state dental boards for additional recommendations and visit: cdc.gov or aami.org

In addition to weekly monitoring, practices should perform biological monitoring:

- ✓ Whenever a new type of packaging material or tray is used
- ✓ After training new sterilization personnel
- ✓ After any sterilizer repair or maintenance
- ✓ After any change in the sterilizer loading process

HELPFUL REMINDERS



Never run the BI on the first load of the day



Always run the BI on a full load if using a table top sterilizer



Sterilant penetration best monitored by placing BIs in the hardest to reach areas of sterilizer

KEY



Place your initials over the date each week that BI testing was performed.



EVERY LOAD: Sure-Check™ Sterilization Pouches

Multi-parameter process indicators will turn from pink to brown indicating Time/Temp/Steam sterilization requirements have been met.



EVERY DAY: SteamPlus™ Integrating Indicator

Use a check mark to indicate daily testing with an integrating indicator. When the dark bar enters the SAFE zone, non-implant loads can be released while waiting for the biological monitoring results.



EVERY WEEK: ConFirm™ Mail-in or In-Office Biological Monitoring System

Indicate Pass (P), Fail (F) or Invalid (I). Also indicate mail-in (M) or in-office (O) biological monitoring type.

JANUARY

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)-(M/O)
				1	2	3	
4	5	6	7	8	9	10	
11	12	13	14	15	16	17	
18	19	20	21	22	23	24	
25	26	27	28	29	30	31	

FEBRUARY

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)-(M/O)
1	2	3	4	5	6	7	
8	9	10	11	12	13	14	
15	16	17	18	19	20	21	
22	23	24	25	26	27	28	

MARCH

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)-(M/O)
1	2	3	4	5	6	7	
8	9	10	11	12	13	14	
15	16	17	18	19	20	21	
22	23	24	25	26	27	28	
29	30	31					

APRIL

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)-(M/O)
			1	2	3	4	
5	6	7	8	9	10	11	
12	13	14	15	16	17	18	
19	20	21	22	23	24	25	
26	27	28	29	30			

MAY

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)-(M/O)
					1	2	
3	4	5	6	7	8	9	
10	11	12	13	14	15	16	
17	18	19	20	21	22	23	
24	25	26	27	28	29	30	
31							

JUNE

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)-(M/O)
	1	2	3	4	5	6	
7	8	9	10	11	12	13	
14	15	16	17	18	19	20	
21	22	23	24	25	26	27	
28	29	30					

JULY

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)-(M/O)
			1	2	3	4	
5	6	7	8	9	10	11	
12	13	14	15	16	17	18	
19	20	21	22	23	24	25	
26	27	28	29	30	31		

AUGUST

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)-(M/O)
						1	
2	3	4	5	6	7	8	
9	10	11	12	13	14	15	
16	17	18	19	20	21	22	
23	24	25	26	27	28	29	
30	31						

SEPTEMBER

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)-(M/O)
		1	2	3	4	5	
6	7	8	9	10	11	12	
13	14	15	16	17	18	19	
20	21	22	23	24	25	26	
27	28	29	30				

OCTOBER

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)-(M/O)
				1	2	3	
4	5	6	7	8	9	10	
11	12	13	14	15	16	17	
18	19	20	21	22	23	24	
25	26	27	28	29	30	31	

NOVEMBER

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)-(M/O)
1	2	3	4	5	6	7	
8	9	10	11	12	13	14	
15	16	17	18	19	20	21	
22	23	24	25	26	27	28	
29	30						

DECEMBER

S	M	T	W	T	F	S	WEEKLY RESULTS (P/F/I)-(M/O)
		1	2	3	4	5	
6	7	8	9	10	11	12	
13	14	15	16	17	18	19	
20	21	22	23	24	25	26	
27	28	29	30	31			



Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | HuFriedyGroup.com

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC its affiliates or related companies, unless otherwise noted. HuFriedyGroup.com

©2025 Hu-Friedy Mfg. Co., LLC. All rights reserved. 2511

